

For Research Use Only

LIMS1 Polyclonal antibody

Catalog Number: 55336-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

55336-1-AP

Size:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004987

GeneID (NCBI):

3987

UNIPROT ID:

P48059

Full Name:

LIM and senescent cell antigen-like domains 1

Calculated MW:

37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB: HeLa cells, A549 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yasmin A Kadry	30254023	J Cell Sci	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of LIMS1 antibody (55336-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-LIMS1 transfected A549 cells.



WB result of LIMS1 antibody (55336-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-LIMS1 transfected HeLa cells.