

For Research Use Only

PRKACG Polyclonal antibody

Catalog Number: 55346-1-AP

1 Publications



Basic Information

Catalog Number:

55346-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002732

GeneID (NCBI):

5568

UNIPROT ID:

P22612

Full Name:

protein kinase, cAMP-dependent,
catalytic, gamma

Calculated MW:

40 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB : mouse testis tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Bingjie Liu	38218754	Hum Cell	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 55346-1-AP (PRKACG antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.