

For Research Use Only

PPP1R3F Polyclonal antibody

Catalog Number:55386-1-AP



Basic Information

Catalog Number:

55386-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_033215

GeneID (NCBI):

89801

UNIPROT ID:

Q6ZSY5

Full Name:

protein phosphatase 1, regulatory
(inhibitor) subunit 3F

Calculated MW:

83 kDa

Observed MW:

83 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB,ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : SMMC-7721 cells,

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 55386-1-AP (PPP1R3F Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.