

For Research Use Only

# TGF Beta Receptor I Polyclonal antibody

Catalog Number: 55391-1-AP

4 Publications

## Basic Information

Catalog Number:

55391-1-AP

Size:

320 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_004612

GeneID (NCBI):

7046

UNIPROT ID:

P36897

Full Name:

transforming growth factor, beta receptor 1

Calculated MW:

56 kDa

Observed MW:

56 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse

Cited Species:

human, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Na Tang	28942050	Exp Parasitol	WB
Jingjing Wu	29077259	J Cell Mol Med	WB
Wanming He	30742067	Oncogene	WB

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

