

For Research Use Only

KLHL9 Polyclonal antibody

Catalog Number: 55394-1-AP

1 Publications



Basic Information

Catalog Number:

55394-1-AP

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_018847

GeneID (NCBI):

55958

UNIPROT ID:

Q9P2J3

Full Name:

kelch-like 9 (Drosophila)

Calculated MW:

69 kDa

Observed MW:

80 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, IP

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB : HeLa cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yao Yao	39708321	Cell Rep	WB,IP,IF

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 55394-1-AP (KLHL9 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.