

For Research Use Only

KLHL9 Polyclonal antibody, PBS Only

Catalog Number: 55394-1-PBS



Basic Information

Catalog Number:

55394-1-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_018847

GeneID (NCBI):

55958

UNIPROT ID:

Q9P2J3

Full Name:

kelch-like 9 (Drosophila)

Calculated MW:

69 kDa

Observed MW:

80 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 55394-1-AP (KLHL9 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 55394-1-PBS in a different storage buffer formulation.