

For Research Use Only

# NOV Polyclonal antibody

Catalog Number: 55449-1-AP

1 Publications



## Basic Information

Catalog Number:

55449-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002514

GeneID (NCBI):

4856

UNIPROT ID:

P48745

Full Name:

nephroblastoma overexpressed gene

Calculated MW:

39 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:50-1:500

## Applications

Tested Applications:

IHC, ELISA

Cited Applications:

IHC

Species Specificity:

human

Cited Species:

human

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

IHC : human nephroblastoma tissue, human prostate cancer tissue

## Background Information

### Notable Publications

| Author         | Pubmed ID | Journal        | Application |
|----------------|-----------|----------------|-------------|
| Xiaojian Huang | 31454155  | J Cell Mol Med | IHC         |

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human nephroblastoma using 55449-1-AP (NOV antibody) at dilution of 1:50 (under 40x lens).