

For Research Use Only

# TRPM6 Polyclonal antibody

Catalog Number: 55455-1-AP

1 Publications



## Basic Information

Catalog Number:

55455-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_017662

GeneID (NCBI):

140803

UNIPROT ID:

Q9BX84

Full Name:

transient receptor potential cation channel, subfamily M, member 6

Calculated MW:

232 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:10-1:100

## Applications

Tested Applications:

IF/ICC, ELISA

Cited Applications:

IHC

Species Specificity:

human

Cited Species:

human

Positive Controls:

IF/ICC : HEK-293 cells,

## Background Information

### Notable Publications

| Author         | Pubmed ID | Journal                   | Application |
|----------------|-----------|---------------------------|-------------|
| Maosheng Cheng | 38619620  | Cancer Immunol Immunother | IHC         |

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of HEK-293 cells using 55455-1-AP (TRPM6 antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG.