

For Research Use Only

SEL1L Polyclonal antibody

Catalog Number: 55459-1-AP



Basic Information

Catalog Number:

55459-1-AP

Size:

390 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005065

GeneID (NCBI):

6400

UNIPROT ID:

Q9UBV2

Full Name:

sel-1 suppressor of lin-12-like (C. elegans)

Calculated MW:

89 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

