

For Research Use Only

# DYNC2H1 Polyclonal antibody

Catalog Number: 55473-1-AP

1 Publications



## Basic Information

Catalog Number:

55473-1-AP

Size:

220 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001080463

GeneID (NCBI):

79659

UNIPROT ID:

Q8NCM8

Full Name:

dynein, cytoplasmic 2, heavy chain 1

Calculated MW:

493 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

IF

Species Specificity:

human, mouse

Cited Species:

human

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Laura Vuolo	30320547	Elife	IF

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

