

For Research Use Only

SARM1 Polyclonal antibody

Catalog Number: 55488-1-AP



Basic Information

Catalog Number:

55488-1-AP

Size:

770 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_015077

GeneID (NCBI):

23098

UNIPROT ID:

Q6SZW1

Full Name:

sterile alpha and TIR motif containing
1

Calculated MW:

79 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:50-1:500

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Positive Controls:

IHC : mouse brain tissue,

**Note-IHC: suggested antigen retrieval with
TE buffer pH 9.0; (*) Alternatively, antigen
retrieval may be performed with citrate
buffer pH 6.0**

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

