

For Research Use Only

GATA1 Monoclonal antibody

Catalog Number: 60011-2-Ig



Basic Information

Catalog Number:

60011-2-Ig

Size:

N/A

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1350

GenBank Accession Number:

BC009797

GeneID (NCBI):

2623

UNIPROT ID:

P15976

Full Name:

GATA binding protein 1 (globin transcription factor 1)

Calculated MW:

43 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

5B1A7

Applications

Tested Applications:

WB,ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data