

# IFNAR2 Monoclonal antibody

Catalog Number: 60013-1-Ig

## Basic Information

Catalog Number:

60013-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgA

Immunogen Catalog Number:

AG0191

GenBank Accession Number:

BC002793

GeneID (NCBI):

3455

UNIPROT ID:

P48551

Full Name:

interferon (alpha, beta and omega)  
receptor 2

Calculated MW:

331 aa, 37 kDa

Observed MW:

58-150 kDa

Purification Method:

Caprylic acid/ammonium sulfate  
precipitation

CloneNo.:

1D4D9

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

