

For Research Use Only

# IFNAR2 Monoclonal antibody

Catalog Number: 60013-2-Ig



## Basic Information

Catalog Number:

60013-2-Ig

Size:

N/A

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG0191

GenBank Accession Number:

BC002793

GeneID (NCBI):

3455

UNIPROT ID:

P48551

Full Name:

interferon (alpha, beta and omega)  
receptor 2

Calculated MW:

331 aa, 37 kDa

Purification Method:

Caprylic acid/ammonium sulfate  
precipitation

CloneNo.:

4G8G1

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data