

For Research Use Only

TIRAP Monoclonal antibody

Catalog Number: 60016-1-Ig



Basic Information

Catalog Number:

60016-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0771

GenBank Accession Number:

BC032474

GeneID (NCBI):

114609

UNIPROT ID:

P58753

Full Name:

toll-interleukin 1 receptor (TIR)
domain containing adaptor protein

Calculated MW:

256 aa, 28 kDa

Observed MW:

35 kDa

Purification Method:

Protein G purification

CloneNo.:

4H5E12

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data