

For Research Use Only

CIRBP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60025-3-PBS



Basic Information

Catalog Number:

60025-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0280

GenBank Accession Number:

BC000403

GeneID (NCBI):

1153

UNIPROT ID:

Q14011

Full Name:

cold inducible RNA binding protein

Calculated MW:

172 aa, 19 kDa

Purification Method:

Protein G purification

CloneNo.:

1D11F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51405-1, CIRBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60025-3-PBS. Detection antibody: 60025-4-PBS.



Cytometric bead array standard curve of MP51405-1, CIRBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60025-3-PBS. Detection antibody: 60025-4-PBS. Standard:Ag16702. Range: 0.098-100 ng/mL.