

For Research Use Only

CIRBP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60025-4-PBS



Basic Information

Catalog Number: 60025-4-PBS	GenBank Accession Number: BC000403	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 1153	CloneNo.: 1D8B9
Isotype: IgG1	UNIPROT ID: Q14011	
Immunogen Catalog Number: AG0280	Full Name: cold inducible RNA binding protein	
	Calculated MW: 172 aa, 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51405-1, CIRBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60025-3-PBS. Detection antibody: 60025-4-PBS.



Cytometric bead array standard curve of MP51405-1, CIRBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60025-3-PBS. Detection antibody: 60025-4-PBS. Standard:Ag16702. Range: 0.098-100 ng/mL.