

For Research Use Only

# Neprilysin/CD10 Monoclonal antibody

Catalog Number: 60034-1-Ig



## Basic Information

Catalog Number:

60034-1-Ig

Size:

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG0427

GenBank Accession Number:

BC101658

GeneID (NCBI):

4311

UNIPROT ID:

P08473

Full Name:

membrane metallo-endorpeptidase

Calculated MW:

750 aa, 86 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

4A10H5

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

MME is also names as neprilysin, CALLA, NEP, SFE, CD10 and belongs to the peptidase M13 family. It is a glycoprotein that is particularly abundant in kidney, where it is present on the brush border of proximal tubules and on glomerular epithelium. CD10 has been considered a useful marker in the diagnosis of renal carcinomas, because of its expression in clear cell and papillary renal cell carcinomas and its absence in chromophobe renal cell carcinomas (PMID:15286660).

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

