

For Research Use Only

Neprilysin/CD10 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60034-5-PBS



Basic Information

Catalog Number:

60034-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG0427

GenBank Accession Number:

BC101658

GeneID (NCBI):

4311

UNIPROT ID:

P08473

Full Name:

membrane metallo-endoropeptidase

Calculated MW:

750 aa, 86 kDa

Purification Method:

Protein A purification

CloneNo.:

1F11A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50348-1, MME,CD10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60034-5-PBS. Detection antibody: 68921-1-PBS. Standard:Ag12506. Range: 3.125-100 ng/mL