

For Research Use Only

# PSMG1 Monoclonal antibody

Catalog Number: 60048-1-Ig



## Basic Information

Catalog Number:

60048-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG0303

GenBank Accession Number:

BC003619

GeneID (NCBI):

8624

UNIPROT ID:

O95456

Full Name:

proteasome (prosome, macropain)  
assembly chaperone 1

Calculated MW:

288 aa, 33 kDa

Observed MW:

44 kDa

Purification Method:

Caprylic acid/ammonium sulfate  
precipitation

CloneNo.:

1B6E1

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

