

For Research Use Only

LAMR1,RPSA Monoclonal antibody

Catalog Number: 60068-1-Ig



Basic Information

Catalog Number:

60068-1-Ig

Size:

810 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG1055

GenBank Accession Number:

BC008867

GeneID (NCBI):

3921

UNIPROT ID:

P08865

Full Name:

ribosomal protein SA

Calculated MW:

33 kDa

Observed MW:

40 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

5H3G8

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

