

For Research Use Only

SEPT2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60075-2-PBS



Basic Information

Catalog Number:

60075-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG1946

GenBank Accession Number:

BC014455

GeneID (NCBI):

4735

UNIPROT ID:

Q15019

Full Name:

septin 2

Calculated MW:

41 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A12B6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50982-1, SEPT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60075-1-PBS. Detection antibody: 60075-2-PBS. Standard:Ag1946. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50982-3, SEPT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60075-3-PBS. Detection antibody: 60075-2-PBS. Standard:Ag1946. Range: 0.098-100 ng/mL.