

For Research Use Only

# SEPT2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60075-3-PBS



## Basic Information

Catalog Number:	60075-3-PBS	GenBank Accession Number:	BC014455	Purification Method:	Protein G purification
Size:	1 mg/ml	GeneID (NCBI):	4735	CloneNo.:	5B3D5
Source:	Mouse	UNIPROT ID:	Q15019		
Isotype:	IgG1	Full Name:	septin 2		
Immunogen Catalog Number:	AG1946	Calculated MW:	41 kDa		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50982-2, SEPT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60075-3-PBS. Detection antibody: 60075-4-PBS. Standard:Ag1946. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50982-3, SEPT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60075-3-PBS. Detection antibody: 60075-2-PBS. Standard:Ag1946. Range: 0.098-100 ng/mL