

For Research Use Only

EEF1D Monoclonal antibody, PBS Only



Catalog Number: 60085-2-PBS

Featured Product

Basic Information

Catalog Number:

60085-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0983

GenBank Accession Number:

BC007847

GeneID (NCBI):

1936

UNIPROT ID:

P29692

Full Name:

eukaryotic translation elongation factor 1 delta (guanine nucleotide exchange protein)

Calculated MW:

31 kDa

Observed MW:

35-40 kDa

Purification Method:

Protein A purification

CloneNo.:

3G4C4

Applications

Tested Applications:

WB, IP, Indirect ELISA, IHC

Species Specificity:

human, mouse, rat

Background Information

EEF1D, also named as EF1D and EF 1 delta, belongs to the EF-1-beta/EF-1-delta family. It is a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. EF-1-beta and EF-1-delta stimulate the exchange of GDP bound to EF-1-alpha to GTP. EEF1D is phosphorylated upon DNA damage, probably by ATM or ATR. The calculated molecular weight of EEF1D is a 31 kDa, but the modified protein is about 35-40 kDa. (PMID: 21936567)

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human pancreas cancer slide using 60085-2-Ig (EEF1D Antibody) at dilution of 1:50. This data was developed using the same antibody clone with 60085-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human pancreas cancer slide using 60085-2-Ig (EEF1D Antibody) at dilution of 1:50. This data was developed using the same antibody clone with 60085-2-PBS in a different storage buffer formulation.



IP result of anti-EEF1D (IP:60085-2-Ig, 5ug; Detection:60085-2-Ig 1:1000) with MCF-7 cells lysate 3200ug. This data was developed using the same antibody clone with 60085-2-PBS in a different storage buffer formulation.