

For Research Use Only

# OPN3 Monoclonal antibody

Catalog Number: 60090-1-Ig



## Basic Information

Catalog Number:

60090-1-Ig

Size:

500 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5075

GenBank Accession Number:

BC036773

GeneID (NCBI):

23596

UNIPROT ID:

Q9H1Y3

Full Name:

opsin 3

Calculated MW:

402 aa, 45 kDa

Observed MW:

45 kDa

Purification Method:

Caprylic Acid/Ammonium Sulfate  
Precipitation

CloneNo.:

1E10D10

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

