

For Research Use Only

# PFN2 Monoclonal antibody

Catalog Number: 60094-1-Ig **2 Publications**



## Basic Information

|                           |            |                           |                |                      |                        |
|---------------------------|------------|---------------------------|----------------|----------------------|------------------------|
| Catalog Number:           | 60094-1-Ig | GenBank Accession Number: | BC018049       | Purification Method: | Protein A purification |
| Concentration:            | N/A        | GeneID (NCBI):            | 5217           | CloneNo.:            | 1D9C4                  |
| Source:                   | Mouse      | UNIPROT ID:               | P35080         |                      |                        |
| Isotype:                  | IgM        | Full Name:                | profilin 2     |                      |                        |
| Immunogen Catalog Number: | AG5043     | Calculated MW:            | 140 aa, 15 kDa |                      |                        |

## Applications

Tested Applications:  
ELISA

Cited Applications:  
WB

Species Specificity:  
human, mouse

Cited Species:  
human

## Background Information

### Notable Publications

| Author      | Pubmed ID | Journal          | Application |
|-------------|-----------|------------------|-------------|
| Min Yu      | 29773903  | Oncogene         | WB          |
| Xiaobin Cui | 29238898  | Mol Cell Biochem | WB          |

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

