

For Research Use Only

TWEAKR Monoclonal antibody

Catalog Number: 60104-1-Ig



Basic Information

Catalog Number:

60104-1-Ig

Size:

N/A

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG2689

GenBank Accession Number:

BC002718

GeneID (NCBI):

51330

UNIPROT ID:

Q9NP84

Full Name:

tumor necrosis factor receptor
superfamily, member 12A

Calculated MW:

129 aa, 14 kDa

Purification Method:

Caprylic Acid/Ammonium Sulfate
Precipitation

CloneNo.:

2E10C10

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human heart tissue were subjected to SDS PAGE followed by western blot with 60104-1-Ig (TWEAKR antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.