

For Research Use Only

FGFR2 Monoclonal antibody

Catalog Number: 60106-1-Ig



Basic Information

Catalog Number:

60106-1-Ig

Size:

2400 μ g/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG3701

GenBank Accession Number:

BC039243

GeneID (NCBI):

2263

UNIPROT ID:

P21802

Full Name:

fibroblast growth factor receptor 2

Calculated MW:

709 aa, 79 kDa

Observed MW:

130 kDa

Purification Method:

Protein A purification

CloneNo.:

5E1D5

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

