

For Research Use Only

# DR6 Monoclonal antibody

Catalog Number: 60107-1-Ig



## Basic Information

Catalog Number:

60107-1-Ig

Size:

600 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1126

GenBank Accession Number:

BC005192

GeneID (NCBI):

27242

UNIPROT ID:

O75509

Full Name:

tumor necrosis factor receptor  
superfamily, member 21

Calculated MW:

15 kDa

Observed MW:

72 kDa

Purification Method:

Caprylic acid/ammonium sulfate  
precipitation

CloneNo.:

4D4F8

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

