

For Research Use Only

Fibromodulin Monoclonal antibody, PBS Only



Catalog Number: 60108-2-PBS

Basic Information

Catalog Number: 60108-2-PBS	GenBank Accession Number: BC035281	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 2331	CloneNo.: 2B11D3
Isotype: IgG1	UNIPROT ID: Q06828	
Immunogen Catalog Number: AG4030	Full Name: fibromodulin	
	Calculated MW: 376 aa, 43 kDa	
	Observed MW: 50-60kDa, 43 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, pig, rabbit

Background Information

Fibromodulin (FMOD), an extracellular matrix protein, belongs to the small leucine-rich proteoglycan family. Fibromodulin is a collagen-binding protein present in many types of connective tissues, e.g. cartilage, tendon, skin, sclera and cornea (PMID: 2531085). It is involved in fibrillogenesis, cell adhesion, and cytokine activity modulation (PMID: 15741214). Fibromodulin may participate in the assembly of the extracellular matrix as it interacts with type I and type II collagen fibrils and inhibits fibrillogenesis in vitro. It may also regulate TGF-beta activities by sequestering TGF-beta into the extracellular matrix. The calculated molecular mass of fibromodulin is 43 kDa, extracellular fibromodulin appears as one single band of 59 kDa due to glycosylation (PMID: 2531085; 22138190).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 60108-2-Ig (FMOD antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60108-2-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 60108-2-Ig (FMOD antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60108-2-PBS in a different storage buffer formulation.