

FXYD6 Monoclonal antibody

Catalog Number: 60111-1-Ig

Basic Information

Catalog Number:

60111-1-Ig

Size:

N/A

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG8538

GenBank Accession Number:

BC018652

GeneID (NCBI):

53826

UNIPROT ID:

Q9H0Q3

Full Name:

FXYD domain containing ion transport
regulator 6

Calculated MW:

95 aa, 11 kDa

Purification Method:

Caprylic Acid/Ammonium Sulfate
Precipitation

CloneNo.:

6H10F4

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

