

For Research Use Only

Hepsin Monoclonal antibody

Catalog Number: 60114-1-Ig



Basic Information

Catalog Number:

60114-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_002151

GeneID (NCBI):

3249

UNIPROT ID:

P05981

Full Name:

hepsin

Calculated MW:

45 kDa

Observed MW:

51 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

2D3A9

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Hepsin, a type II membrane-associated protein, is a cell surface proteases. It acts directly by degrading the extracellular matrix of surrounding cells and also functions as an activator of other extracellular proteases. It was reported to be involved in multiple biological processes.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

