

For Research Use Only

CLCA2 Monoclonal antibody

Catalog Number: 60118-2-Ig



Basic Information

Catalog Number:

60118-2-Ig

Size:

N/A

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG5914

GenBank Accession Number:

BC041096

GeneID (NCBI):

9635

UNIPROT ID:

Q9UQC9

Full Name:

chloride channel accessory 2

Calculated MW:

104 kDa

Purification Method:

Caprylic Acid/Ammonium Sulfate
Precipitation

CloneNo.:

4E11B8

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data