

For Research Use Only

# Kv1.2 Monoclonal antibody

Catalog Number: 60120-1-Ig



## Basic Information

Catalog Number:

60120-1-Ig

Size:

500 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5477

GenBank Accession Number:

BC043564

GeneID (NCBI):

3737

UNIPROT ID:

P16389

Full Name:

potassium voltage-gated channel,  
shaker-related subfamily, member 2

Calculated MW:

499 aa, 57 kDa

Observed MW:

75 kDa

Purification Method:

Caprylic Acid/Ammonium Sulfate  
Precipitation

CloneNo.:

1D2D11

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data