

For Research Use Only

# IL-9R Monoclonal antibody

Catalog Number: 60132-1-Ig



## Basic Information

Catalog Number:

60132-1-Ig

Size:

500 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG11033

GenBank Accession Number:

BC051337

GeneID (NCBI):

3581

UNIPROT ID:

Q01113

Full Name:

interleukin 9 receptor

Calculated MW:

522 aa, 57 kDa

Observed MW:

55 kDa, 57 kDa

Purification Method:

Caprylic Acid/Ammonium Sulfate  
Precipitation

CloneNo.:

2A12H5

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data