

For Research Use Only

# CHGA Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60135-3-PBS



## Basic Information

Catalog Number:

60135-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0807

GenBank Accession Number:

BC006459

GeneID (NCBI):

1113

UNIPROT ID:

P10645

Full Name:

chromogranin A (parathyroid secretory protein 1)

Calculated MW:

51 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3H11B7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50241-1, CHGA /Chromogranin A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 688S1-1-PBS. Detection antibody: 60135-3-PBS. Standard:Ag19780. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50241-2, CHGA/Chromogranin A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60135-2-PBS. Detection antibody: 60135-3-PBS. Standard:Ag19780. Range: 0.781-100 ng/mL