

For Research Use Only

# GSTO1 Monoclonal antibody

Catalog Number: 60150-1-Ig



## Basic Information

Catalog Number: 60150-1-Ig	GenBank Accession Number: BC000127	Purification Method: Protein A purification
Size: N/A	GeneID (NCBI): 9446	CloneNo.: 1E3E9
Source: Mouse	UNIPROT ID: P78417	
Isotype: IgM	Full Name: glutathione S-transferase omega 1	
Immunogen Catalog Number: AG7199	Calculated MW: 28 kDa	

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

