

For Research Use Only

GABRD Monoclonal antibody

Catalog Number: 60163-1-Ig



Basic Information

Catalog Number: 60163-1-Ig	GenBank Accession Number: BC033801	Purification Method: Protein G purification
Size: 500 µg/ml	GeneID (NCBI): 2563	CloneNo.: 4E9D9
Source: Mouse	UNIPROT ID: O14764	
Isotype: IgG1	Full Name: gamma-aminobutyric acid (GABA) A receptor, delta	
Immunogen Catalog Number: AG8064	Calculated MW: 51 kDa	

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human

Background Information

Gamma-aminobutyric acid (GABA) is the major inhibitory neurotransmitter in the mammalian brain. The GABRD gene encodes the delta subunit of the ligand-gated chloride channel for gamma-aminobutyric acid (GABA). This protein migrates as a 55-65kd protein in SDS-PAGE.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data