

GABRD Monoclonal antibody

Catalog Number: 60163-1-Ig

Basic Information

Catalog Number:

60163-1-Ig

Size:

500 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8064

GenBank Accession Number:

BC033801

GeneID (NCBI):

2563

UNIPROT ID:

O14764

Full Name:

gamma-aminobutyric acid (GABA) A
receptor, delta

Calculated MW:

51 kDa

Purification Method:

Protein G purification

CloneNo.:

4E9D9

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Background Information

Gamma-aminobutyric acid (GABA) is the major inhibitory neurotransmitter in the mammalian brain. The GABRD gene encodes the delta subunit of the ligand-gated chloride channel for gamma-aminobutyric acid (GABA). This protein migrates as a 55-65kd protein in SDS-PAGE.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

