

For Research Use Only

KCNH7 Monoclonal antibody

Catalog Number: 60166-1-Ig



Basic Information

Catalog Number:

60166-1-Ig

Size:

500 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4723

GenBank Accession Number:

BC035815

GeneID (NCBI):

90134

UNIPROT ID:

Q9NS40

Full Name:

potassium voltage-gated channel,
subfamily H (eag-related), member 7

Calculated MW:

1196 aa, 135 kDa

Observed MW:

120-130 kDa

Purification Method:

Protein G purification

CloneNo.:

2F2G1

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

