

For Research Use Only

CD122/IL2RB Monoclonal antibody

Catalog Number: 60174-2-Ig



Basic Information

Catalog Number: 60174-2-Ig	GenBank Accession Number: BC025691	Purification Method: Protein G purification
Size: N/A	GeneID (NCBI): 3560	CloneNo.: 5G1E6
Source: Mouse	ENSEMBL Gene ID: ENSG00000100385	Recommended Dilutions: WB 1:500-1:1000
Isotype: IgG2a	UNIPROT ID: P14784	
Immunogen Catalog Number: AG3929	Full Name: interleukin 2 receptor, beta	
	Calculated MW: 551 aa, 61 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Recombinant protein protein,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 60174-2-Ig (IL2RB antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Recombinant protein were subjected to SDS PAGE followed by western blot with 60174-2-Ig (IL2RB antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.