

For Research Use Only

# CD122/IL2RB Monoclonal antibody

Catalog Number: 60174-2-Ig



## Basic Information

Catalog Number:

60174-2-Ig

Size:

N/A

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG3929

GenBank Accession Number:

BC025691

GeneID (NCBI):

3560

ENSEMBL Gene ID:

ENSG00000100385

UNIPROT ID:

P14784

Full Name:

interleukin 2 receptor, beta

Calculated MW:

551 aa, 61 kDa

Purification Method:

Protein G purification

CloneNo.:

5G1E6

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Recombinant protein protein,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 60174-2-Ig (IL2RB antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Recombinant protein were subjected to SDS PAGE followed by western blot with 60174-2-Ig (IL2RB antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.