

For Research Use Only

# CD122/IL-2RB Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60174-3-PBS



## Basic Information

Catalog Number:

60174-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3929

GenBank Accession Number:

BC025691

GeneID (NCBI):

3560

ENSEMBL Gene ID:

ENSG00000100385

UNIPROT ID:

P14784

Full Name:

interleukin 2 receptor, beta

Calculated MW:

551 aa, 61 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3F12D9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50451-1, CD122/IL2RB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60174-3-PBS. Detection antibody: 60174-1-PBS. Standard:Ag3929. Range: 0.098-25 ng/mL.



Cytometric bead array standard curve of MP50451-3, CD122/IL2RB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60174-2-PBS. Detection antibody: 60174-3-PBS. Standard:Ag3929. Range: 0.098-3.125 ng/mL.