

For Research Use Only

CD3E Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60181-2-PBS

Basic Information

Catalog Number: 60181-2-PBS	GenBank Accession Number: BC049847	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 916	CloneNo.: 5B12E11
Source: Mouse	ENSEMBL Gene ID: ENSG00000198851	
Isotype: IgG2b	UNIPROT ID: P07766	
Immunogen Catalog Number: AG11797	Full Name: CD3e molecule, epsilon (CD3-TCR complex)	
	Calculated MW: 207 aa, 23 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50222-1, CD3E Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60181-2-PBS. Detection antibody: 60181-1-PBS. Standard:Ag11797. Range: 0.098-50 ng/mL



Cytometric bead array standard curve of MP50222-2, CD3E Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60181-2-PBS. Detection antibody: 60181-3-PBS. Standard:Ag11797. Range: 0.098-100 ng/mL