

For Research Use Only

CD3E Monoclonal antibody, PBS Only (Detector)



Catalog Number: 60181-3-PBS

Basic Information

Catalog Number:

60181-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG11797

GenBank Accession Number:

BC049847

GeneID (NCBI):

916

ENSEMBL Gene ID:

ENSG00000198851

UNIPROT ID:

P07766

Full Name:

CD3e molecule, epsilon (CD3-TCR complex)

Calculated MW:

207 aa, 23 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3C12G12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50222-2, CD3E Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60181-2-PBS. Detection antibody: 60181-3-PBS. Standard:Ag11797. Range: 0.098-100 ng/mL.