

For Research Use Only

SMAD4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60182-2-PBS

Basic Information

Catalog Number:

60182-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0299

GenBank Accession Number:

BC002379

GeneID (NCBI):

4089

UNIPROT ID:

Q13485

Full Name:

SMAD family member 4

Calculated MW:

60 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C2H7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50149-1, SMAD4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60182-2-PBS. Detection antibody: 60182-3-PBS. Standard:Ag0299. Range: 0.098-50 ng/mL