

For Research Use Only

GFAP Monoclonal antibody, PBS Only

Catalog Number: 60190-2-PBS



Basic Information

Catalog Number:

60190-2-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG10452

GenBank Accession Number:

BC013596

GeneID (NCBI):

2670

UNIPROT ID:

P14136

Full Name:

glial fibrillary acidic protein

Calculated MW:

432 aa, 50 kDa

Purification Method:

Protein A purification

CloneNo.:

3C4D10

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

