

For Research Use Only

GFAP Monoclonal antibody

Catalog Number: 60190-3-Ig **1 Publications**



Basic Information

Catalog Number: 60190-3-Ig	GenBank Accession Number: BC013596	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 2670	CloneNo.: 3B2D11
Isotype: IgG1	UNIPROT ID: P14136	
Immunogen Catalog Number: AG10452	Full Name: glial fibrillary acidic protein	
	Calculated MW: 432 aa, 50 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IF

Species Specificity:
human

Cited Species:
mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yi-Fu Zhu	39021082	Zool Res	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data