

For Research Use Only

GFAP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60190-4-PBS



Basic Information

Catalog Number:

60190-4-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10452

GenBank Accession Number:

BC013596

GeneID (NCBI):

2670

UNIPROT ID:

P14136

Full Name:

glial fibrillary acidic protein

Calculated MW:

432 aa, 50 kDa

Purification Method:

Protein G purification

CloneNo.:

5E2F6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP80011-1, GFAP Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 81063-1-PBS. Detection antibody: 60190-4-PBS.



Cytometric bead array standard curve of MP80011-1, GFAP Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 81063-1-PBS. Detection antibody: 60190-4-PBS. Standard:Ag10452. Range: 0.391-400 ng/mL.