

For Research Use Only

# DNAJB2 Monoclonal antibody

Catalog Number: 60204-1-Ig **1 Publications**



## Basic Information

Catalog Number: 60204-1-Ig	GenBank Accession Number: BC011609	Purification Method: Protein A purification
Concentration: 770 µg/ml	GeneID (NCBI): 3300	CloneNo.: 2G8C5
Source: Mouse	UNIPROT ID: P25686	Recommended Dilutions: WB: 1:500-1:2000
Isotype: IgG2a	Full Name: DnaJ (Hsp40) homolog, subfamily B, member 2	
Immunogen Catalog Number: AG1271	Calculated MW: 36 kDa	
	Observed MW: 80 kDa	

## Applications

Tested Applications: WB, ELISA	Positive Controls: WB : SH-SY5Y cells,
Cited Applications: WB	
Species Specificity: human, mouse	
Cited Species: mouse	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yuxin Qing	40035989	Inflammation	WB

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 60204-1-Ig (DNAJB2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.